



# Sculptor

## Anti-aging System



TECHNIQUE CREATES BEAUTY

ALL YOU NEED, ALL IN ONE



## Sculptor

In the facial rejuvenation treatment, we need to target the skin multi-level comprehensive treatment, Sculptor is mainly on the epidermis to the fascia layer, accelerate the epidermal metabolism, stimulate the dermis in the contraction of elastic fibers and collagen newborn, and tighten the fascia, to achieve the effect of tightening wrinkles.



MORE THAN  
15 YEARS



PRESENT IN  
115  
COUNTRIES



SCIENTIFICALLY-P  
ROVEN  
155 SCIENTIFIC

# SCULPTOR

## Anti-aging system

## Leading Technology


Traditional bipolar radiofrequency energy is transmitted from the positive electrode to the negative electrode and penetrates the skin shallowly.

Traditional monopolar radiofrequency energy is transmitted from the positive electrode to the negative electrode plate. Due to the uneven distribution of impedance in the human body, the energy can only pass through the superficial layer of the skin.

Based on many years of radio frequency technology and clinical experience, we successfully launched this epoch-making product. Unlike previous unipolar RF, Sculptor does not require a negative plate to be attached to the human body, making it more convenient to use. Radiofrequency energy is emitted downwards in a vortex with the treatment head as the center, and can penetrate deeply into all layers of the skin. Sculptor can automatically match the impedance of human skin to achieve a comfortable thermal effect from the inside out.


Using high-frequency 13MHZ Sculptor technology, the electromagnetic field causes the quantum particles to move with each other (invisible to the naked eye), causing thermal shock, accelerating mutual friction to generate heat, and penetrating into the subcutaneous tissue in a vortex-like manner. Because its heat is diffused, it can be heated. All layers of the skin achieve completely non-invasive, all-round three-dimensional anti-aging.

Accelerate metabolism and blood circulation, stimulate fibroblasts to secrete collagen and elastin, make the skin resistant, restore healthy skin, improve skin texture and color, make it rosy, shiny, firm and elastic.




# The development history of physical beauty and non-invasive anti-aging.

Technology paints a magnificent poem of physical beauty



## 2023 KES 13.56MHZ Sculptor





# A Full Spectrum of Solutions

- ◆ Face lifting
- ◆ Skin lifting around eyes
- ◆ Skin lifting on neck
- ◆ Improve conditions of sagging skin
- ◆ Forehead wrinkle removal
- ◆ Degerming and acne removal
- ◆ Cellulite reduction and slimming
- ◆ skin tightening and skin whitening

## Clinic Advantages

- 1.Layered cell-molecule friction heating technology penetrates deeply and achieves three-dimensional anti-aging.
- 2.It only heats the molecular gaps and does not damage the inherent form of the substance, which is safer.
- 3.Emits energy in the form of a vortex, diffusing energy from the inside out, transferring the energy evenly without burning the skin.
- 4.High comfort, no anesthesia required during treatment, no wounds.
- 5.The heat stays for a long time, and the heat will stay in the skin for a period of time after the operation is completed, continuing to stimulate the regeneration of collagen.
- 6.The large hand can operate on the face, the small hand can operate around the eyes, and the negative pressure hand can better operate on looser skin or neck.
- 7.It lasts for a long time and the effect is remarkable. One course of treatment can last for 1-1.5 years.

## Before and After



# Clinic Advantages

- 13.56MHZ RF Technology
- The maximum power can reach 120W
- No negative electrode required
- Automatic impedance matching
- Comfortable and no tingling sensation
- 5 Levels temperature Adjustable



# Handle Advantage

- Screen display function (gear, countdown, treatment head)
- Gears adjustable Circularly
- Changeable treatment tip Automatic Identification
- Patented hand plug design
- Vacuum technology lifting function
- LED reminder light
- The treatment head adopts advanced insulation coating



# Advantages of independent research and development of appearance

- Novel appearance
- 13.3-inch screen (1920x1080 resolution)
- Background setting function (energy, negative pressure adjustable)
- Automatic identification and power-on self-test of treatment Tips
- Countdown treatment time
- OEM customized cart and appearance acceptable



# Technique Specification

Machine Consumption	150 VA (Max)
Power Supply	220V /110V 50Hz/60Hz
Display screen size	13.3 inch
Screen Resolution	1920x1080
RF Frequency	13Mhz
MAX Output power	120W
HANDPIECE TIPS	I (eye around , check, forehead) II (FACE , neck, Jaw)
Energy	1-5 Levels Adjustable
Setting	RF energy、 Vacuum intensity

## About us

Beijing KES Biology Technology Co., Ltd is one of global leading suppliers of suppliers of Diode Laser, IPL, Velashape, Nd-yag Q Switch Laser, CO2 Fractional Laser, Cryolipolysis solutions for the aesthetic and medical markets.



## Factory

Tel : +86-10- 84405668

Email : [info@keslaser.com](mailto:info@keslaser.com)

Web : [www.keslaser.com](http://www.keslaser.com)

 No. 15, Section 2, Mafang Industrial Park, Pinggu Park, Pinggu District, Beijing 101200 China